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
Datasheet

Cios Select

Introducing flat-detector capabilities to routine surgery

siemens.com/cios-select

SIEMENS
Healthineers 

Table of contents

	Page
System specification	
C-arm	3
Collimator system	3
PC Hardware	3
X-ray generator/tube	
2.3 kW High-frequency generator	4
Single tank with stationary-anode tube	4
Flat detector 21 cm x 21 cm (8.3" x 8.3")	5
X-ray release switches	5
Handswitch	5
Footswitch	5
Monitor cart	6
Clinical workflow	
Patient data administration	7
Exam preparation	7
Image acquisition	7
Operating modes	7
Fluoroscopy	7
Single Image	7
Subtraction ¹⁾ /Roadmap ¹⁾	7
Intelligent Dose Efficiency Algorithm	7
CARE program (Combined Applications to Reduce Exposure)	
CAREVISION	8
CAREPROFILE	8
Dose optimization	8

	Page
Image display/processing	
Image display	9
Image processing	9
Digital Density Optimization (DDO)	9
Text/graphic functions	9
Data transfer and documentation	
DICOM network interfaces	10
WLAN ¹⁾	11
CD/DVD	11
Printer interface ¹⁾	11
USB export	11
DVI interface for external monitor ¹⁾	11
Administration and service	
Administration	11
Service ¹⁾	11
Room planning	
Operating data	12
Environmental conditions (operation)	12
Dimensions and weight	12
Dimensions in cm (inches)	13
Legal notice	14

¹⁾ Option

Table of contents

	Page		Page
System specification		Image display/processing	
C-arm	3	Image display	9
Collimator system	3	Image processing	9
PC Hardware	3	Digital Density Optimization (DDO)	9
X-ray generator/tube		Text/graphic functions	9
2.3 kW High-frequency generator	4	Data transfer and documentation	
Single tank with stationary-anode tube	4	DICOM network interfaces	10
Flat detector 21 cm x 21 cm (8.3" x 8.3")	5	WLAN ¹⁾	11
X-ray release switches	5	CD/DVD	11
Handswitch	5	Printer interface ¹⁾	11
Footswitch	5	USB export	11
Monitor cart	6	DVI interface for external monitor ¹⁾	11
		Administration and service	
Clinical workflow		Administration	11
Patient data administration	7	Service ¹⁾	11
Exam preparation	7		
Image acquisition	7	Room planning	
Operating modes	7	Operating data	12
Fluoroscopy	7	Environmental conditions (operation)	12
Single Image	7	Dimensions and weight	12
Subtraction ¹⁾ /Roadmap ¹⁾	7	Dimensions in cm (inches)	13
Intelligent Dose Efficiency Algorithm	7	Legal notice	14
CARE program			
(Combined Applications to Reduce Exposure)			
CAREVISION	8		
CAREPROFILE	8		
Dose optimization	8		

System specification

Page 4 – 7

Clinical workflow

Page 8 – 11

Accessories

Page 12

Room planning

Page 13 – 15

C-arm

C-arm with mechanical brakes

2.8.1

Orbital movement

130 ° (– 40° to + 90°)

Angulation

± 190°

Horizontal movement

20 cm (7.9")

Immersion depth

73 cm (28.7")

Swivel range

± 12°

Vertical travel

42 cm (16.5"), motorized

Source to image-receptor distance

100 cm (39.4")

Free space

78 cm (30.7"), with attached grid

Collimator system

Iris diaphragm (lead)

For concentric, radiation-free collimation

Semitransparent slot diaphragm (lead)

For symmetric, radiation-free collimation, with unlimited rotation

PC Hardware

Image acquisition system

Intel-compatible dual-core microprocessor with PCI bus architecture, Windows 7, 64 bit, 4 GB RAM, S-ATA drive, USB 2.0, high-performance and interface cards for the detector/X-ray system

Acquisition memory on hard disk

Storage and postprocessing of all acquired images in a common patient folder

Hard-disk capacity

150,000 images on hard disk irrespective of matrix size

Power supply

2 TB

Integrated, uninterruptible power supply helps ensure that image and patient data are secure in the event of a power outage

* Option

System specifications

X-ray generator/tube

2.3 kW High-frequency generator

Nominal peak output power	2.3 kW
Inverter control frequency	29 kHz to 44 kHz
kV range	40 kV to 110 kV
Single Image	3 mA to 24 mA
Continuous fluoroscopy	0.2 to 13 mA
Pulsed fluoroscopy	3 mA to 24 mA
	Pulse width: 7 ms to 40 ms
	Pulse rate: Variable frame rate 0.5 to 15 f/s 30 f/s in continuous fluoro mode
Power and dose management	Easy selection of suitable power and dose levels for each clinical application

Single tank with stationary-anode tube

Focal spot nominal value	0.6 / 1.0
Nominal voltage	110 kV
Anode heat dissipation	40,800 J/min 55,000 HU/min
Anode heat storage capacity (IEC 60613)	75,000 J 101,000 HU
Optical anode angle	9°
Inherent filtration (IEC 60601)	3 mm Al with 70 kVp / 0.1 mm Cu
Single-tank heat dissipation	3,400 J/min 4,594 HU/min
Single-tank heat storage capacity (physical)	800,000 J 1,100,000 HU
Continuous heat dissipation	80 W
Max. uninterrupted fluoro time	50 min at 300 W

¹ Option

System specifications

Flat detector 21 cm x 21 cm (8.3" x 8.3")

Amorphous silicon flat detector based on indirect conversion technology

High-performance digital imaging system

Input fields (active field)	Mag 0 (full format) 20.58 cm x 20.58 cm (8.1" x 8.1")	Mag 1 15 cm x 15 cm (5.9" x 5.9")	Mag 2 10 cm x 10 cm (3.9" x 3.9")
Material	aSi with CsI scintillator		
Pixel size	205 µm		
Matrix	1004 x 1004 pixels		
Digitization depth	16 bits		
Detective Quantum Efficiency DQE(0), typical, RQAS	80 % at 0 lp/mm 65 % at 1 lp/mm 40 % at 2 lp/mm 22 % at 2.43 lp/mm (Nyquist frequency)		
Modulation Transfer Function (MTF), typical	55 % at 1 lp/mm 22 % at 2 lp/mm		
Resolution on monitor with 21 x 21 FD (measured in accordance with DIN 6868-150)	Overview (Mag 0) – acquisition	2.43 lp/mm	
	Overview (Mag 0) – fluoroscopy	2.43 lp/mm	
	Format switchover (Mag 1) – acquisition	2.43 lp/mm	
	Format switchover (Mag 1) – fluoroscopy	2.43 lp/mm	
	Format switchover (Mag 2) – acquisition	2.43 lp/mm	
	Format switchover (Mag 2) – fluoroscopy	2.43 lp/mm	
Anti-scatter grid (detachable)	Pb 17 : 1, 70 lines/cm, $f_0 = 100$ cm		
X-ray release switches			
Handswitch	Handswitch with 3 m spiral cable, 1 x X-ray release button, 1 x image store button		
Footswitch	Standard footswitch for radiation release (5 m cable)		
	Multifunctional footswitch wireless ¹⁾ with advanced functionality		
	Transmission: Bluetooth BT4.0LT Standard (ISM-Band 2,402 GHz – 2,480 GHz, frequency hopping).		
	Battery type 3 x LR14 (Type C, Baby) Duration: min. 517 days with 1 h per day accumulated use and 8000 mAh battery capacity Footswitch IP-class IPx8		

¹⁾ Option

China 3E-61

Monitor cart

Displays

19" TFT high-contrast, high-luminance color display

Screen diagonal	19" (48 cm)
Image display	1280 x 1024 (pixels)
Brightness (DICOM calibrated), typical	400 cd/m ²
Maximum brightness	800 cd/m ²
Horizontal / vertical viewing angle	176 ° / 176 °
Contrast ratio, typical	1000 : 1
Backlight LED technology	

19" TFT color display

Screen diagonal	19" (48 cm)
Image display	1280 x 1024 (pixels)
Brightness, typical	180 cd/m ²
Maximum brightness	300 cd/m ²
Horizontal / vertical viewing angle	178 ° / 178 °
Contrast ratio, typical	2000 : 1

Interface description

Video splitter*	2 single-link DVI-D; female connectors; up to 1920 x 1200 @ 60 Hz or 1080 p
USB	Data interface USB 2.0 Charging interface USB 2.0
LAN	Gigabit
Power cable	7 m
Monitor cart connection cable	8.5 m

1-6-1

Clinical workflow

Patient data administration

Patient registration

Retrieval of patient list and examination data from the hospital/radiology information system (HIS/RIS)
Emergency patient registration
Study and image data administration
Configurable patient registration

Exam preparation

Applications Manager

Dedicated, application-related user programs
Administration and selection of applications and application groups

Image acquisition

Operating modes

Selection of application-specific fluoroscopy and radiography curves for the individual operating modes

Fluoroscopy

Image acquisition 0.5 f/s to 15 f/s; 30 f/s with continuous fluoro mode
Image storage: all images, every n^{th} image
Digital filtration
Moving weighted averaging for low-noise image display with minimum dose
Optimized gray-scale visualization based on image analysis

Single Image

Digital filtration
Optimized gray-scale visualization based on image analysis

Subtraction ¹⁾/Roadmap ¹⁾

Image acquisition 0.5 f/s to 15 f/s; 30 f/s with continuous fluoro mode
Image storage: all images, every n^{th} image
Subtraction angiography
Peak Opacification for iodine contrast (MaxOp) and CO₂ contrast (MinOp)
Anatomical landmarking from 0 % to 100%
Manual pixelshift function to correct subtraction runs
Filtering of mask and fill images for contrast enhancement
Simultaneous dual-channel output for image acquisition and postprocessing, simultaneous storage of fill image
Roadmap technique, to position a catheter precisely in a blood vessel under fluoroscopy

Intelligent Dose Efficiency Algorithm

The Intelligent Dose Efficiency Algorithm is a smart dose management system specially designed for the Cios C-arm family
The result is an always optimized balance of image quality and dose as well as automatic contrast and brightness regulation

¹ Option

Clinical workflow

CARE program (Combined Applications to Reduce Exposure)

CARE is a Siemens Healthineers initiative to reduce radiation dose. Cios Select is equipped with state-of-the-art features to reduce radiation dose to patients, surgeons and staff including automatic dose management. The result is an always optimized balance of image quality and dose.

CAREVISION

Variety of dedicated Exam Sets to adjust the settings for the right dose, e.g. for pediatric cases (low) or obese patients (high)
Pulsed fluoroscopy with a pulse rate of up to 15 p/s
Easy selection of dose levels and operating modes including dedicated low-dose programs

CAREPROFILE

Radiation-free positioning of primary collimators through graphical display in the LIH image on the image monitor

Dose optimization

Integrated dose measuring chamber with automatic transfer of the accumulated dose into a radiation report
Selection of dose levels
Additional copper filter for further reduction of the patient dose
System-integrated FD Laser aimer^{*)}
Detachable grid, e.g. for pediatric applications

^{*)} Option

Clinical workflow

Data transfer and documentation

DICOM network interfaces

DICOM Send/Storage Commitment ¹⁾	DICOM interface for image data communication in a clinical network (PACS) based on the DICOM 3 standard Sending, receiving and storing of images Archiving confirmation from the image archive
DICOM Print ¹⁾	For printing within the network, on a DICOM-compatible camera or DICOM-compatible printer
DICOM Query/Retrieve ¹⁾	Retrieval of studies from a digital archive, a workstation, or other imaging systems; e.g. MR, CT Multi-modality viewing
DICOM Worklist/MPPS ¹⁾	Get Worklist function for importing patient data from a data management system (RIS/HIS). XA, CR and DX worklist entries supported, configurable Modality Performed Procedure Step (MPPS) function for sending examination statistics and dose information to a data management system
DICOM Dose Structured Report ¹⁾	Sending of dose values for each study to an archiving system
DICOM Advanced ¹⁾	DICOM Advanced contains all the functions of DICOM Dose Structured Report, plus: DICOM Send/Storage Commitment DICOM Print DICOM Query/Retrieve DICOM Worklist/MPPS

¹⁾ Option

Clinical workflow

Image display/processing

Image display

Aspect ratio 5:4, corresponding to 1280 x 1024 matrix, 1k x 1k image content
 Split screen (1, 16 on 1)
 Digital zoom, fixed zoom, roaming
 Magnification (detector zoom)
 Digital image rotation
 Movie function for playback of scenes and auto replay function
 Digital shutters
 Horizontal and vertical image reversal
 LSH (Last Scene Hold)

Image processing

Application-specific lookup tables (LUTs) for optimum contrast and brightness
 Spatial frequency filtration for edge-enhanced image display
 Edge enhancement
 Noise reduction
 Motion detection with active noise reduction
 Metal correction

Digital Density Optimization (DDO)

Digital Density Optimization reduces the dynamic range of an image, allowing the contrast of structures to be emphasized without loss of information in bright and dark image areas

Text/graphic functions

Text: annotation, image comments, R/L marking
 Graphics: quantification with distance and angle measurements

Digital video imaging system

High-resolution X-ray video system in maintenance-free CCD technology with 1024 x 1024 (1K²) full-size sensor
 Automatic gain control for constant image brightness

High contrast and high spatial resolution

1.3.4

Video matrix

1K²

Digital image rotation

± 360° continuous

1.3.6

Clinical workflow

Data transfer and documentation

WLAN ¹⁾

WLAN client module for wireless transmission of DICOM image data, e.g. to a PACS

Compatible with 802.11 a/b/e/g/h/i/n WLAN standards

Operation within the 2.4 / 5 GHz frequency bands

802.11 i, 802.1 x, WPA/WPA2

WPA2 Enterprise supplicants EAP-TLS, EAP-TTLS, (MSCHAPv2), EAPP-EAP (MSCHAPv2) as security/authentication features

Supports TKIP and AES for data encryption

Supports DHCP client

CD/DVD

DVD drive for digital image storage on CD-R, DVD+R or DVD-R for offline data exchange in DICOM 3, TIFF and AVI formats

Printer interface ¹⁾

Digital printers for printing on paper or paper/film

USB export

For digital image storage to a USB device in DICOM, TIFF and AVI formats

DVI interface for external monitor ¹⁾

Live monitor (A):

Video splitter output for connecting an external live monitor

Reference monitor (B):

Video splitter output for connecting an external reference monitor

DVI interface (splitter), without galvanic isolation

Administration and service

Administration

HIPAA ¹⁾

X-ray passcode ¹⁾

Service ¹⁾

Remote Desktop (TeamViewer)

Room planning

Operating data

Power requirements	100 V, 110 V, 120 V, 127 V, 200 V, 230 V, 240 V, ($\pm 10\%$), 50/60 Hz (± 1 Hz)
Unit fuse protection (internal)	100 V to 127 V 20 A slow-blow fuse 200 V to 240 V 15 A slow-blow fuse
Maximum power consumption	3.8 KW
Standby power consumption	713 W
Voltage / Current values	Continuous 14 A (110 V) / 7 A (230 V) Short-time 25 A (110 V) / 16 A (230 V)
Internal line impedance	Ri max. 0.3 ohms for 100 V to 127 V Ri max. 0.8 ohms for 200 V to 240 V

Environmental conditions (operation)

Temperature range	+ 10 °C to + 35 °C
Relative humidity	20 % to 75 %, non-condensing
Barometric pressure	700 hPa to 1060 hPa

Dimensions and weight


Chassis (L x w x h)	182.5 cm x 80 cm x 159.8 cm (71.9" x 31.5" x 62.9")
Monitor cart (L x w x h)	68.6 cm x 73.6 cm x 179.1 cm (27" x 29.1" x 70.5")

Weight of C-arm chassis without options

Cios Select	275 kg (605 lbs)
Monitor cart	150 kg (330 lbs)

Room planning

Dimensions in cm (inches)

On account of certain regional limitations of sales rights and service availability, we cannot guarantee that all products/services/features included in this brochure are available through the Siemens sales organization worldwide. Availability and packaging may vary by country and are subject to change without prior notice. Some/All of the features and products described herein may not be available in the United States or other countries.

The information in this document contains general technical descriptions of specifications and options as well as standard and optional features that do not always have to be present in individual cases.

Siemens reserves the right to modify the design, packaging, specifications and options described herein without prior notice. Please contact your local Siemens sales representative for the most current information.

Cios Select is currently under development; is not for sale in the U.S. Its future availability cannot be guaranteed.

In the interest of complying with legal requirements concerning the environmental compatibility of our products (protection of natural resources and waste conservation), we may recycle certain components where legally permissible.

Using the same extensive quality assurance measures as for factory-new components, we guarantee the quality of these recycled components.

Note: Any technical data contained in this document may vary within defined tolerances. Original images always lose a certain amount of detail when reproduced.

Caution: Federal law restricts this device to sale by or on the order of a physician.

Cios Select 14455694
VA20

Siemens Healthineers Headquarters
Siemens Healthcare GmbH
Henkestr. 127
91052 Erlangen
Germany
Phone +49 9131 84-0
siemens.com/healthineers

Legal Manufacturer
Siemens Shanghai
Medical Equipment Ltd.
No. 278, Zhouzhu Road
Shanghai, 201318
P. R. China